**Vakuum und Destillationscontroller DC1 / DC1P**  
**Vacuum and Distillation Controller DC1 / DC1P**  
**Appareil de contrôle du vide et de la distillation DC1 / DC1P**

**Bedienungsanleitung**  
**Operating Instructions**  
**Mode d'emploi**

<b>1</b>	<b>Allgemeine Hinweise .....</b>	<b>3</b>	<b>1</b>	<b>General Information.....</b>	<b>3</b>	<b>1</b>	<b>Généralités .....</b>	<b>3</b>
<b>2</b>	<b>Beschreibung.....</b>	<b>4</b>	<b>2</b>	<b>Description.....</b>	<b>4</b>	<b>2</b>	<b>Description .....</b>	<b>4</b>
<b>3</b>	<b>Installation.....</b>	<b>5</b>	<b>3</b>	<b>Installation .....</b>	<b>5</b>	<b>3</b>	<b>Installation .....</b>	<b>5</b>
3.1	Anschlüsse.....	5	3.1	Connectors.....	5	3.1	Connecteurs .....	5
3.1.1	Steuerausgänge .....	5	3.1.1	Control Outputs .....	5	3.1.1	Sorties de commande .....	5
3.1.2	RS232 (opt.) .....	5	3.1.2	RS232 (opt.).....	5	3.1.2	RS232 (opt.) .....	5
3.2	Inbetriebnahme.....	6	3.2	Initial Start-Up .....	6	3.2	Première Mise en Service.....	6
<b>4</b>	<b>Bedienung .....</b>	<b>6</b>	<b>4</b>	<b>Operation .....</b>	<b>6</b>	<b>4</b>	<b>Fonctionnement .....</b>	<b>6</b>
4.1	Einschalten.....	6	4.1	Power On .....	6	4.1	Mise sous tension .....	6
4.2	Konfiguration .....	7	4.2	Configuration.....	7	4.2	Configuration .....	7
4.3	Manual Mode ( nur DC1/1 ) .....	8	4.3	Manual Mode ( nur DC1/1 ) .....	8	4.3	Mode Manuel (DC1/1).....	8
4.4	Programm Mode ( nur DC1P ) .....	10	4.4	Program Mode ( nur DC1P ).....	10	4.4	Mode Programmation (DC1P) .....	10
4.5	Auto Mode .....	13	4.5	Auto Mode.....	13	4.5	Mode automatique .....	13
4.6	Tempopär Mode .....	15	4.6	Tempopary Mode .....	15	4.6	Mode Temporaire .....	15
4.7	Serielle Schnittstelle .....	16	4.7	Serial Interface .....	16	4.7	Interface Série .....	16
<b>5</b>	<b>Wartung und Pflege .....</b>	<b>16</b>	<b>5</b>	<b>Maintenance.....</b>	<b>16</b>	<b>5</b>	<b>Maintenance .....</b>	<b>16</b>
<b>6</b>	<b>Zubehör .....</b>	<b>16</b>	<b>6</b>	<b>Accessories .....</b>	<b>16</b>	<b>6</b>	<b>Accessoires.....</b>	<b>16</b>
<b>7</b>	<b>Technische Daten .....</b>	<b>17</b>	<b>7</b>	<b>Technical Data .....</b>	<b>17</b>	<b>7</b>	<b>Données techniques .....</b>	<b>17</b>

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## 1 General Information

The operation manual has been written for the DC1/1 and DC1P vacuum and distillation controller.

This instrument represents latest technical standards and is manufactured according to established safety regulations. Nevertheless, improper use leads to measurement errors and damage of the device.

Therefore it is mandatory to read and follow this operating instructions.

### Delivery:

The following items are delivered:

- DC1 vacuum controller with integrated aerate valve
- 2 hose nozzles with o-ring sealing
- operating instructions, certificate of conformity
- counter-plugs for control outputs

### Application:

The DC1 vacuum and distillation controller with internal ceramics sensor is suitable for measuring and controlling pressure, e.g. in evaporators or central vacuum supplies.

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

- 1 Display Operation Mode
- 2 Display Setpoint
- 3 Display Actual Pressure
- 4 Power Switch
- 5 SET-key
- 6 Arrow-keys
- 7 Start-keys for automatic and manual mode
- 8 STOP key
- 9 Socket for power supply
- 10 Control outputs, 2x female, 6pole  
for vacuum valve, pump or cooling water
- 11 Aeration Valve
- 12 Mounting Bolt, with M4 thread
- 13 Hose Nozzle 8mm, vacuum
- 14 Valve Block
- 15 Hose Nozzle 8mm, aeration

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


#### 3.1 Connectors

##### 3.1.1 Control Outputs

Maximum switching power for each  
valve: 5 W !

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

### 3.2 Initial Start-Up

Connect cable for mains power supply and vacuum hose nozzle at the back of the instrument to the process chamber. Connect the necessary control outputs for vacuum valve (chapter 3.1.1) or frequency converter for pump speed control and optionally the switching outputs for pump or cooling water (chapter 3.1.2).

## 4 Operation

### 4.1 Power On

The DC1 can be operated in different modes in order to enable optimum controlling for various applications.

The power switch is located at the right side of the instrument.

After switching on the DC1 the currently set pressure unit (hPa=mbar or Torr) is displayed in the right window named "Actual Pressure". Afterwards the corresponding displays show the current operation mode, the current setpoint and the actual pressure.

Controlling in the respective operation mode ist started with key "AUTO" resp. "MAN".

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## 4.2 Configuration

To change basic settings of the DC1 in the configuration mode,

**AUT** hold the "AUTO"-key pressed when switching the instrument on.

In the display "Actual Pressure" the current display unit appears. It can be switched between "hPa" (mbar) or "Torr" using the arrow keys.

**SET** press Set-key,

in the Mode-display "A1" is shown, right of it the numeric value of the automatic-parameter, e.g. 38 (see chapter 4.5). By means of the arrow keys it can be set from 20 (sensitive) to 50 (insensitive).

**SET** press Set-key,

in the Mode-display "Pb" is shown, right of it the numeric value of the P-control-parameter Xp: 1% (soft) to 50% (hard).

**STOP** With the "STOP"-key you can leave the configuration mode at any time.

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

### 4.3 Manual Mode (DC1/1 only)

In the Manual Mode a pressure setpoint is controlled, which is set manually. You can adjust the following setpoint-profile:

- 1: quick decrease to setpoint
- 2: ramp-profile, gradient of decrease according to pressure decrease and decrease time
- 3: process end, valves closed

The following parameters can be adjusted:

- setpoint, 1-999
- hysteresis (h), 1-100
- pressure decrease (d), 0-100
- decrease time (t), toff (without) - 99

The dialog display for the input is the "Setpoint"-Display. Parameters can be set using the arrow keys. After 10s without any keypress, the instrument switches back to display the actual pressure.

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## Input example:

- SET** setpoint is displayed, 250 (hPa)
- SET** hysteresis is displayed, h\_5 (hPa)
- SET** pressure decrease is displayed, d\_45 (hPa)
- SET** decrease time is displayed, t\_15 (min)
- SET** actual pressure is displayed.

## Start and stop process:

- MAN** start controlling, the Mode-display shows "H1" (LED "MAN" is on)

If a pressure decrease is programmed, the remaining decrease time is displayed alternating with the operation mode "H1". When the process is finished the Mode-display shows "End" and all valves are closed.

If the "MAN"-key is pressed a second time during the process, the process chamber is evacuated with maximum pump speed.

- STOP** stop controlling, aerate process chamber.

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
#### 4.4 Program Mode ( DC1P only )

With model DC1P you can program set-point-profiles. Therefore 9 programs with max. 9 steps/each can be edited and stored.

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
Parameter values can be adjusted by using the arrow-keys. If there no input is done, programming is stopped after 10s. Programming also ends when the "STOP"-key is pressed or when finally the parameter value for hysteresis (HYS) is entered.

In an arbitrary program step the ramp time or hold time can be set to "tend", which defines this step as program end.

Input example:

<b>SET</b>	"Pn"/"1", selection of program number
<b>SET</b>	"A1"/"1", ramp time step 1
<b>SET</b>	"SP1"/"500", setpoint step 1
<b>SET</b>	"t1"/"2", hold time step 1
	...
<b>SET</b>	"A6"/"tend", ramp time step 6 program end
<b>SET</b>	"HYS"/"5", hysteresis

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Start and stop program control:

**MAN** start program control, (LED "MAN" is on)

The Mode-display shows program- and current step-number, e.g. "P1.2" (prog. 1, step 2). During the hold phases the program/step number is displayed alternating with the remaining hold time.

If the "MAN"-key is pressed a second time during the program run, the process chamber is evacuated with maximum pump speed.

During program control only the active program can be edited.

At the end of the program the Mode-display shows "End" and all valves are closed.

**STOP** stop program, aerate process chamber.

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#### 4.5 Auto Mode

In the Automatic Mode the first boiling point ("bP": boiling point) is found automatically during pressure decrease. For optimization you can adjust the automatic parameter (see chapter 4.1). A small numeric value achieves a more sensitive automatic. If the parameter is set too small, a premature detection of the boiling point can be the result. Setting the parameter too high can cause overriding of the boiling point.

With model DC1/1 you can set a pressure decrease ramp (see chapter 4.2) - just like in the Manual-mode. This controlled decrease is done after the first boiling point has been detected.

Hereby the following parameters are adjustable:

- setpoint (of no importance)
- hysteresis (h), 1-100
- pressure decrease (d), 0-100
- decrease time (t), toff (without) - 99

Start and stop process:

**AUT** start controlling, the Mode-display shows "A1", LED AUTO is on.

Arrow- and SET-keys are locked. When the first boiling point is detected, "bP" appears in the Mode-display.

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If -with model DC1/1- a pressure decrease is programmed, the remaining decrease time is displayed alternating with "bP". When the process is finished the Mode-display shows "End" and all valves are closed.

With model DC1P the detected boiling point is controlled, until you stop the process pressing the STOP-key.

**STOP** stop controlling,  
aerate process chamber.

After detection of the boiling point you can switch to Manual-mode at any time. To do this:

**MAN** press "MAN"-key, the instrument switches into Manual-mode.

**SET** By pressing the "SET"-key the boiling point "bP" is stored as new setpoint.

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#### 4.6 Temporary Mode

In the modes "Manual" and "Program" the user is able to intervene at any time independent of the programmed process controlling by setting a temporary setpoint, which is not stored as operating parameter:

By a short keypress on an arrow-key the actual pressure is set as

temporary setpoint (LED "MAN" is flashing).

The temporary setpoint can now be adjusted using the arrow-keys.

If you keep pressing the key, the process chamber is evacuated with maximum pump speed.

Hereby the setpoint is steadily adapted to the current actual pressure.

If you keep pressing the key, the aeration valve is opened shortly 3-times, afterwards keeps opened.

Hereby the setpoint is steadily adapted to the current actual pressure.

With the "MAN"-key you can leave the temporary mode.

The original setpoint is controlled again. Ending the process with the "STOP"-key is also possible at any time.

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## 4.7 Serial Interface

Controllers of type DC1 are optionally equipped with a serial interface RS232.

If the 9-pole Sub-D-socket is connected to a PC by the interface cable (accessory), measurement data can be recorded and visualized via PC (VacuGraph Windows-software, see accessories).

### Interface Parameters:

2400 Baud, 1 startbit, 8 data bits, 1 stop bit, no parity, no handshake;

Communication corresponds to Thyracont-interface-protocol, transmitting ASCII-code.

## 5 Maintenance

Please handle with care. The front can be cleaned by a little bit wet cloth.

## 6 Accessories

The following accessories are available:

- External vacuum valve, nickel coated, 1,5m cable (PTFE valve on request)
- Stand holder
- Software VacuGraph for graphic data display and measurement via PC.

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**7 Technical Data**

Range	1400 - 1 mbar 1050 - 1 Torr
Resolution	1 mbar
Control Output	Relais 230 VAC, 1 A
Materials with Media-Contact	ceramics, PP, Viton
Measuring Rate	20 ms
Operating Temperature	0 ... 40°C
Storage Temperature	-10 ... +60°C
Power Supply	95 - 265 VAC 50/60 Hz
Power Consumption	max. 13 W
Analog Output (optional)	0 - 10 VDC potential free max. 10 kOhm
Control Outputs	24 VDC, max. 5W
Connection f. Vacuum/Aeration	G1/8 or hose nozzle
Display	LED 13mm
Dimensions	145x85x50(85) mm
Weight	550g

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**Thyracont Electronic GmbH**