Vakuum und Destillationscontroller DC1S
Vacuum and Distillation Controller DC1S
Appareil de contrôle du vide et de la distillation DC1S

Bedienungsanleitung
Operating Instructions
Mode d'emploi

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1 General Information

The operation manual has been written for the DC1S digital vacuum and distillation controller.

This instrument represents latest technical standards and is manufactured according to established safety regulations. Nevertheless, improper use lead to measurement errors and damage of the device.

Therefore it is mandatory to read and follow this operating instructions.

Delivery:

The following items are delivered:

- DC1S vacuum controller with integrated aerate valve
- 3 hose nozzles with o-ring sealing
- operating instructions, certificate of conformity
- counter-plugs for control outputs

Application:

The DC1S vacuum and distillation controller with internal ceramic sensor and continuous valves for vacuum and venting is suitable for measuring and controlling pressure, e.g. in evaporators or central vacuum supplies.

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

- 1 Display Operation Mode
- 2 Display Setpoint
- 3 Display Actual Pressure
- 4 SET-key
- 5 Arrow-keys
- 6 Start-keys for automatic and manual mode
- 7 Power Switch
- 8 STOP key
- 9 Port for venting
- 10 Port f. vacuum chamber, G1/4 + hose nozzle 8/12mm
- 11 Valve block
- 12 Port for vacuum, G1/4 + hose nozzle 8/12mm
- 13 Mains socket
- 14 Control outputs, 2x female, 6pole for pump or cooling water
- 15 RS232

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





3.2 Initial Start-Up

Connect the cable for mains supply and your vacuum source, process chamber and –if needed- the pipe for gas inlet to the corresponding ports on the rear side of the DC1S.

Connect optionally the switching outputs for pump or cooling water (chapter 3.1.2).

4 Operation

4.1 Power On

The DC1S can be operated in different modes in order to enable optimum controlling for various applications.

The power switch is located at the bottom right side of the front panel.

After switching on the DC1S the currently set pressure unit (hPa=mbar or Torr) is displayed in the right window named "Actual Pressure". Afterwards the corresponding displays show the current operation mode, the current setpoint and the actual pressure.

Controlling in the respective operation mode is started with key "AUTO" resp. "MAN".

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4.2 Configuration

To change basic settings of the DC1S in the configuration mode,

AUT

hold the "AUTO"-key pressed when switching the instrument on.

In the display "Actual Pressure" the current display unit appears. It can be switched between "hPa" (mbar) or "Torr" using the arrow keys.

SET

press Set-key,

in the Mode-display "A1" is shown, right of it the numeric value of the automatic-parameter, e.g. 38 (see chapter 4.5). By means of the arrow keys it can be set from 20 (sensitive) to 50 (insensitive).

STOP

With the "STOP"-key you can leave the configuration mode at any time.

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

4.3 Manual Mode

In the Manual Mode a pressure setpoint is controlled, which is set manually. You can adjust the following setpoint-profile:

- 1: quick decrease to setpoint
- 2: ramp-profile, gradient of decrease according to pressure decrease and decrease time
- 3: process end, valves closed

The following parameters can be adjusted:

- setpoint, 1-999
- pressure decrease (d), 0-100
- decrease time (t), toff (without) - 99

The dialog display for the input is the "Setpoint"-Display. Parameters can be set using the arrow keys. After 10s without any keypress, the instrument switches back to display the actual pressure.

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Input example:

- SET** setpoint is displayed, 250 (hPa)
- SET** pressure decrease is displayed, d_45 (hPa)
- SET** decrease time is displayed, t_15 (min)
- SET** actual pressure is displayed.

Start and stop process:

- MAN** start controlling, the Mode-display shows "H1" (LED "MAN" is on)

If a pressure decrease is programmed, the remaining decrease time is displayed alternating with the operation mode "H1". When the process is finished the Mode-display shows "End". After 3 more Minutes the process chamber is vented automatically.

If the "MAN"-key is pressed a second time during the process, the process chamber is evacuated with maximum pump speed.

- STOP** stop controlling, aerate process chamber.

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4.5 Auto Mode

In the Automatic Mode the first boiling point ("bP": boiling point) is found automatically during pressure decrease. For optimization you can adjust the automatic parameter (see chapter 4.1). A small numeric value achieves a more sensitive automatic. If the parameter is set too small, a premature detection of the boiling point can be the result. Setting the parameter too high can cause overriding of the boiling point.

With model DC1S you can set a pressure decrease ramp (see chapter 4.2) - just like in the Manual-mode. This controlled decrease is done after the first boiling point has been detected.

Hereby the following parameters are adjustable:

- setpoint (of no importance)
- pressure decrease (d), 0-100
- decrease time (t), toff (without) - 99

Start and stop process:

AUT

start controlling, the Mode-display shows "A1", LED AUTO is on.

Arrow- and SET-keys are locked. When the first boiling point is detected, "bP" appears in the Mode-display.

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If -with model DC1S- a pressure decrease is programmed, the remaining decrease time is displayed alternating with "bP". When the process is finished the Mode-display shows "End". After 3 more Minutes the process chamber is vented automatically.

With model DC1SP the detected boiling point is controlled, until you stop the process pressing the STOP-key.

STOP stop controlling,
aerate process chamber.

After detection of the boiling point you can switch to Manual-mode at any time. To do this:

MAN press "MAN"-key, the instrument switches into Manual-mode.

SET By pressing the "SET"-key the boiling point "bP" is stored as new setpoint.

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4.7 Serial Interface

Controllers of type DC1S are optionally equipped with a serial interface RS232. If the 9-pole Sub-D-socket is connected to a PC by the interface cable (accessory), measurement data can be recorded and visualized via PC (VacuGraph Windows-software, see accessories).

Interface Parameters:

2400 Baud, 1 startbit, 8 data bits, 1 stop bit, no parity, no handshake;
Communication corresponds to Thyracont-interface-protocol, transmitting ASCII-code.

5 Maintenance

Please handle with care. The front can be cleaned by a little bit wet cloth.

6 Accessories

The following accessories are available:

- Stand holder
- Software VacuGraph for graphic data display and measurement via PC.

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7 Technical Data

Range	1400 - 1 mbar 1050 - 1 Torr
Max. Overload	max. 2bar abs.
Resolution	1 mbar
Accuracy	0,5% f.s.
Materials with Media-Contact	ceramics, PP, Viton
Measuring Rate	20 ms
Operating Temperature	0 ... 40°C
Storage Temperature	-10 ... +60°C
Power Supply	95 - 265 VAC 50/60 Hz
Power Consumption	max. 13 W
Analog Output (optional)	0 - 10 VDC potential free max. 10 kOhm
Control Output	24 VDC, max. 5W
Connection f. Vacuum/Aeration	G1/4 or hose nozzle
Display	LED 13mm
Dimensions	155x143x115 mm
Weight	1100g

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Declaration of Conformity

This declaration is valid for following product:

Equipment: Vacuum controller

Type: DC1 S

Hereby the equipment is confirmed to comply with the requirements set out in the Council Directive on the Approximation of the Laws of the Member States relating to Electromagnetic Compatibility (89/336/EEC) and the Council Directive relating to Low Voltage 73/23/EEC.

The following company is responsible for this declaration:

Thyracont Electronic GmbH
Max-Emanuel-Strabe 10
94036 Passau

The measurements were carried out in accredited laboratories.

For the evaluation of above mentioned Council Directives for Electromagnetic Compatibility and for Low Voltage following standards were consulted:

DIN EN61326: 2002-03

DIN EN 55022: 2003-09

DIN EN 61000-4-2: 2001-12, DIN EN 61000-4-3: 2001-12, DIN EN 61000-4-4: 2002-07,
 DIN EN 61000-4-5: 2001-12, DIN EN 61000-4-6: 2001-12,
 DIN EN 61000-4-8: 2001-12, DIN EN 61000-4-11: 2001-12

DIN EN 61000-3-2: 2001-12
 DIN EN 61000-3-3: 2002-05

DIN EN 61010: 2002-08

22.09.04, Passau

Date / Place

Signature of responsible Person